



Protocols

1. DNA extraction
2. PCR
3. Sequence analysis

1. DNA extraction

Material is shock-frozen in liquid nitrogen, and stored at -80°C . The frozen material is ground using a high-throughput disruptor (TissueLyser, Qiagen, Germany), and the DNA is extracted using the InnuPrep Plant DNA Kit (Macherey-Nagel) following the instructions of the producer. Around 50-100 mg of plant material is sufficient. The kit is based on affinity to a column, whereby the DNA binds to the matrix, such that contaminants can be eluted. Subsequently, the DNA can be eluted by elution buffer or nuclease-free water. Typical volumes are 100 μl of eluate. Quality is controlled by spectrophotometry via the Nano-Drop. The device is first initialised by pipetting 1.5 μl of nuclease-free water. Then the device is calibrated by pipetting 1.5 μl of the blank (the solvent for the DNA), then the DNA sample is measured (1.5 μl). The values give the concentration of DNA (in $\text{ng}/\mu\text{l}$) and its purity (determined as A_{260}/A_{280} measuring contamination by proteins and A_{260}/A_{230} measuring contamination by polysaccharides). The A_{260}/A_{280} should be ideally 1.8, higher values indicate RNA contamination, lower values indicate protein contamination. The A_{260}/A_{230} value should be between 1.8 and 2.2. The concentration should then be adjusted to 50 $\text{ng}/\mu\text{l}$ for optimal use in genomic PCR.

2. PCR

Sample volume is 10 μl . Each sample with 1 μl buffer (10x), 1 μl BSA (10 mg/ml), 0.2 μl of each primer (10 pmol), 0.2 μl dNTPs (10 mmol each), and 0.1 μl Taq (New England Biolabs, Frankfurt) and ~ 0.5 μl template (50 $\text{ng}/\mu\text{l}$).

Table 1: PCR protocol for amplification of *tnrH-psbA*

step		temperature	Time
initial denaturation		95°C	2:00 min
denaturation	35x	95°C	0:30 min
annealing		60°C	0:30 min
elongation		68°C	1:00 min
final elongation		68°C	5:00 min
storage		12°C	-

Table 2: oligonucleotide primers for typical plant barcodes.

target	Sequence	Amplicon (bp)
rbcl ^a *	F: 5-ATGTCACCACAAACAGAGACTAAAGC-3 R: 5-CGTGGTGGACTTGATTTTAC-3	599
matK	F: 5-ACCCAGTCCATCTGGAAATCTTGGTTC-3 R: 5-CGTACAGTACTTTTGTGTTTACGAG-3	878
trnL-F igs	F: 5-CGAAATCGGTAGACGCTACG-3 R: 5-ATTTGAACTGGTGACAGAG-3	803-1209
trnH-psbA igs**	F: 5-GTTATGCATGAACGTAATGCTC-3 R: 5-CGCGCATGGTGGATTACAATCC-3	420-446

*source for the oligonucleotide primers for rbcl^a

Levin RA, Wagner WL, Hoch PC, et al. (2003) Family Level Relationships of Onagraceae Based on Chloroplast rbcl and ndhF Data. *Am J Bot* 90, 107-115

**source for the oligonucleotide primers for trnH-psbA

Sang T et al. (1997) Chloroplast DNA phylogeny, reticulate evolution, and biogeography of *Paeonia* (Paeoniaceae). *Am J Bot* 84, 1120.

Tate JA, Simpson BB (2003): Paraphyly of *Tarasa* (Malvaceae) and diverse origins of the polyploid species. *Systematic Botany* 28, 723–737.

Trouble shooting for recalcitrant samples (e.g. dried material from commercial samples). Here, the DNA can be extracted using the Invisorb® MSpin Plant Mini Kit (Strattec Molecular) following the instructions of the producer. Further, for the PCR 10 mg/ml of bovine serum albumin were added.

3. Sequence analysis

The sequences are processed like explained in the sequence tutorial. Basically, the quality of the sequences is assessed by visual inspection of the chromatograms, by BLAST search of the sequence and validation of the sequence identity, and by alignment of the two reads (forward and reverse complement) using the Alignment Explorer of the freeware MEGA (<http://www.megasoftware.net/>). In addition, the NCBI database can be searched by BLAST for homologues, which can be selected by help of the Taxonomy View tool to get a relevant collection of homologues from related taxa. The alignment can then be trimmed and used to construct and bootstrap a Neighbor-Joining (NJ) tree.