



## Protocols

1. Microscopic diagnostics
2. DNA extraction
3. PCR
4. Sequence analysis

## 1. Microscopic diagnostics

---

If conventional fresh samples are difficult to be investigated, one can discolorate the tissue by chloral hydrate maceration:

**NOTE: Chloral hydrate is toxic and irritant. You have to work under a fume hood, wear gloves and protective glasses**

- Hand sections on a slide.
- Drain water, add 3 droplets of chloral hydrate, add coverslip
- Heat up specimen three times, let it cool
- If there are air bubbles caught, drain in chloral hydrate from the side (without heating)
- Carefully drain excess fluid with filter paper

Images can be viewed under 20x magnification at the Leica microscope in 516 (M420, Leica; Bensheim) using a digital camera (DFC 500, Leica; Bensheim). The area of the visual field is 439  $\mu\text{m}$  x 439  $\mu\text{m}$ .

---

## 2. DNA extraction

---

Material is shock-frozen in liquid nitrogen, and stored at  $-80^{\circ}\text{C}$ . The frozen material is ground using a high-throughput disruptor (TissueLyser, Qiagen, Germany), and the DNA is extracted using the InnuPrep Plant DNA Kit (Macherey-Nagel) following the instructions of the producer. Around 50-100 mg of plant material is sufficient. The kit is based on affinity to a column, whereby the DNA binds to the matrix, such that contaminants can be eluted. Subsequently, the DNA can be eluted by elution buffer or nuclease-free water. Typical volumes are 100  $\mu\text{l}$  of eluate. Quality is controlled by spectrophotometry via the Nano-Drop. The device is first initialised by pipetting 1.5  $\mu\text{l}$  of nuclease-free water. Then the device is calibrated by pipetting 1.5  $\mu\text{l}$  of the blank (the solvent for the DNA), then the DNA sample is measured (1.5  $\mu\text{l}$ ). The values give the concentration of DNA (in  $\text{ng}/\mu\text{l}$ ) and its purity (determined as  $A_{260}/A_{280}$  measuring contamination by proteins and  $A_{260}/A_{230}$  measuring contamination by polysaccharides). The  $A_{260}/A_{280}$  should be ideally 1.8, higher values indicate RNA contamination, lower values indicate protein contamination. The  $A_{260}/A_{230}$  value should be between 1.8 and 2.2. The concentration should then be adjusted to 50  $\text{ng}/\mu\text{l}$  for optimal use in genomic PCR.

### 3. PCR

Sample volume is 10 µl. Each sample with 1 µl buffer (10x), 1 µl BSA (10 mg / ml), 0.2 µl of each primer (10 pmol), 0.2 µl dNTPs (10 mmol each), and 0.1 µl Taq (New England Biolabs, Frankfurt) and ~0.5 µl template (50 ng/µl).

**Table 1:** PCR protocol for amplification of trnH-psbA

step		temperature	Time
initial denaturation		95°C	2:00 min
denaturation	35x	95°C	0:30 min
annealing		60°C	0:30 min
elongation		68°C	1:00 min
final elongation		68°C	5:00 min
storage		12°C	-

**Table 2:** oligonucleotide primers for typical plant barcodes.

target	Sequence	Amplicon (bp)
rbcLa*	F: 5-ATGTCACCACAAACAGAGACTAAAGC-3 R: 5-CGTGGTGGACTTGATTTTAC-3	599
matK	F: 5-ACCCAGTCCATCTGGAAATCTTGGTTC-3 R: 5-CGTACAGTACTTTTGTGTTTACGAG-3	878
trnL-F igs	F: 5-CGAAATCGGTAGACGCTACG-3 R: 5-ATTTGAACTGGTGACAGAG-3	803-1209
trnH-psbA igs**	F: 5-GTTATGCATGAACGTAATGCTC-3 R: 5-CGCGCATGGTGGATTCCACAATCC-3	420-446

\*source for the oligonucleotide primers for rbcLa

Levin RA, Wagner WL, Hoch PC, et al. (2003) Family Level Relationships of Onagraceae Based on Chloroplast rbcL and ndhF Data. *Am J Bot* 90, 107-115

\*\*source for the oligonucleotide primers for trnH-psbA

Sang T et al. (1997) Chloroplast DNA phylogeny, reticulate evolution, and biogeography of *Paeonia* (Paeoniaceae). *Am J Bot* 84, 1120.

Tate JA, Simpson BB (2003): Paraphyly of *Tarasa* (Malvaceae) and diverse origins of the polyploid species. *Systematic Botany* 28, 723–737.

**Trouble shooting for recalcitrant samples (e.g. dried material from commercial samples).** Here, the DNA can be extracted using the Invisorb® MSpin Plant Mini Kit (Stratec Molecular) following the instructions of the producer. Further, for the PCR 10 mg/ml of bovine serum albumin were added.

---

#### 4. Sequence analysis

---

The sequences are processed like explained in the sequence tutorial. Basically, the quality of the sequences is assessed by visual inspection of the chromatograms, by BLAST search of the sequence and validation of the sequence identity, and by alignment of the two reads (forward and reverse complement) using the Alignment Explorer of the freeware MEGA (<http://www.megasoftware.net/>). In addition, the NCBI database can be searched by BLAST for homologues, which can be selected by help of the Taxonomy View tool to get a relevant collection of homologues from related taxa. The alignment can then be trimmed and used to construct and bootstrap a Neighbor-Joining (NJ) tree.