

CPS page ruler, Thermo Fisher

pre stained marker

Gel type	Tris-Glycine				Tris-Acetate		Bis-Tris					
Gel concentration	4-20%	8-16%	10-20%	8%	10%	12%	15%	3-8%	7%	4-12%	10%	12%
Running buffer	Tris-Glycine				Tris-Acetate		Bis-Tris					
	Apparent Molecular Weights, kDa											
10			130		170	170	170		100			
20	170	170	130	170	130	100	100		100			
30	150	130	100	120	100	70	70		120	140	140	140
40	100	100	70	100	70	55	55	150	100	115	115	115
50	70	70	55	70	55	40	40	120	85	80	80	80
60	55	55	40	55	40	35	35	85	65	65	65	65
70	40	40	35	40	35	25	25	65	50	50	50	50
80	35	35	25	35	25	15	15	50	40	40	40	40
90	25	25	15	25	15	10	10	40	30	30	30	30
100	15	15	10	15	10	10	10	30	25	25	25	25
	10	10	10	10	10	10	10	25	25	10	10	10

